HOW A **NEUTRON BEAM FACILITY WORKS**

Overview of an experiment

Researchers come from all over Canada and the world to probe materials with neutron beams to find solutions to challenges in health, industry, and science. From the initial inquires into the feasibility of an experiment to the final interpretation of results, CNBC scientists and technical staff provide support to our users to ensure that this national resource is accessible to any user.

As neutrons pass through a material, the material changes the properties of the beam, such as the direction, energy, or magnetic polarization. By detecting these changes, researchers can deduce certain properties of the material such as atomic structure or stress.



Data analysis typically continues after a user travels back to their home institution. CNBC scientists follow up with the users to assist with the analysis and interpretation of results.



Sample Material

A sample material is placed in the emerging

beam. As the neutrons pass through, the material changes the properties of the

beam, such as the direction, energy, and

scattered in many directions. A chamber around the material controls conditions such

magnetic polarization. Typically, the beam is

as temperature, pressure, or magnetic fields.

User Arrival

CNBC, typically after preparing a sample of a cases, samples may be prepared on site.

A researcher travels to the material for study. In some

Several tubes through the reactor wall allow some neutrons to exit in the shape of a beam. Excess neutrons are absorbed by the reactor wall.

Beam Production

Beam Preparation

A crystalline material diffracts the beam, that is, it divides the beam according to the energies of the neutrons. A channel is positioned to allow only neutrons of a certain desired energy to proceed to the sample material. The remaining neutrons are absorbed in the wall of the large cylinder encasing the crystal.

Neutron Production Neutrons are tiny particles

that reside in atoms. When

uranium atoms are split in

the reactor core, neutrons

are released in every


direction with a large

spectrum of energies.



Neutron Detection

A mobile detector system determines the intensity of the scattered beam in various directions.







User Interface and Electronics

A specialized electronic system controls each portion of the beam line and collects the experimental data. Work stations provide the user interface to control the experiment and perform preliminary data analysis.